







# WNS – Labels

Document revision: 01

|                     | Name           | Function             | Date       | Signature                                                                            |
|---------------------|----------------|----------------------|------------|--------------------------------------------------------------------------------------|
| Author              | Elisa Veronese | Inventis Quality Eng | 2022.03.03 |   |
| Verified & Approved | Diego Fiorin   | Inventis R&D Eng     | 2022.03.03 |  |

## WNS Labels

The device does not include any variants. Nonetheless it can be stuck with different labels according to the medical device it is sold with: three different Synapsys medical devices can be connected to Synapsys WNS, i.e. Synapsys VNG, Synapsys DVNS, and Synapsys Nystalab. For traceability purposes, the WNS label allows understanding if the camera has been sold as part of a VNG (2 in the following table), as part of a DVNS (3), as part of a Nystalab medical devices (4), or as spare part (5). The labels (1, 6, 7) do not change according to the medical device.

|   | Internal labels                                                                     | Description                                                                                                                                      |
|---|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 |   | <p>Always present on WNS (Dimensions: 37.5x44 mm).</p> <p>It contains the <b>manufacturer's information</b> (i.e. name and contact address).</p> |
| 2 |  | <p>Used when the WNS is sold as part of a Synapsy VNG medical device (Dimensions: 30x12.5 mm).</p>                                               |
| 3 |  | <p>Used when the WNS is sold as part of a Synapsy D-VNS medical device (Dimensions: 30x12.5 mm).</p>                                             |
| 4 |  | <p>Used when the WNS is sold as part of a Synapsy NystaLab medical device (Dimensions: 30x12.5 mm).</p>                                          |
| 5 |  | <p>Used when the WNS is sold as a spare part (Dimensions: 30x12.5 mm).</p>                                                                       |

|   |                                                                                                                                                                                                                                                                                                                                                                                          |                                                                       |
|---|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|
| 6 | <div style="border: 1px solid black; padding: 5px; width: fit-content; margin: auto;"> <p>This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.</p> </div> | As per FCC requirements. Always present on WNS (Dimensions: 15x17 mm) |
| 7 |                                                                                                                                                                                                                                                                                                         | Always present on WNS.                                                |



View from the top, battery inserted. Protection tip for USB cable removed (left) and inserted (right)



View from the top, battery and protection tip for USB cable removed.



View from the patient's side: detail on the label in the bottom part




View from the patient's side: detail on the label in the upper part.



View from the bottom.


# Radio module label

|   | Internal labels                                                                                                                                      | Description                                  |
|---|------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------|
| 1 | <div><p><b>SYNAPSYS SAS</b><br/>2 rue Marc Donadille<br/>13013 Marseille - France</p><p>FCC ID: 2A4C2-WNSXX</p><p>M/N:<br/>CC3135MODRNMMOB</p></div> | Manufacturer information on the radio module |



Label applied on the module (left). Board embedding the module mounted (the module cannot be shown, being on the opposite part of the board)

# External Labels

|   | External labels                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | Description                                                                  |
|---|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------|
| 1 | <div><p><b>Made in France</b></p><p><b>TRANSPORT AND STORAGE:</b></p><p> -10°C — +55°C</p><p> 10% — 95%</p><p> 700hPa — 1060hPa</p><p> <b>SYNAPSYS SAS - 2 rue Marc Donadille - 13013 Marseille (France)</b></p></div> | Manufacturer information and Transport and storage environmental conditions. |

|   |                                                                                   |    |    |    |    |         |    |    |
|---|-----------------------------------------------------------------------------------|----|----|----|----|---------|----|----|
| 2 |  | AT | BE | BG | CH | CY      | CZ | DE |
|   |                                                                                   | DK | EE | EL | ES | FI      | FR | HR |
|   |                                                                                   | HU | IE | IT | IS | LI      | LT | LU |
|   |                                                                                   | LV | MT | NL | NO | PL      | PT | RO |
|   |                                                                                   | SE | SI | SK | TR | UK (NI) |    |    |

As required by European RED directive.

## External Labels Applied on the Packaging

