






Senstream Research Ring Application
E-label exhibit
01/30/23

Home Page



The screenshot shows the 'Record Session' screen of the Senstream Research Ring Application. At the top, the status bar displays the time 5:07, signal strength, Wi-Fi, and battery icons. The app header includes the time and a profile icon on the left, and an 'Edit' link on the right. The main title 'Record Session' is prominently displayed. Below it, the 'RINGS' section contains a button labeled 'Add Senstream Ring'. The 'SESSION DETAILS' section features a 'Session Type' field with a dropdown menu showing 'Scoring: PPG 200Hz, EDA 5Hz, Temperature 5Hz'. Below this are input fields for 'Session Name' and 'Notes'. A large button at the bottom of the details section reads 'Add a Senstream Ring to Continue'. The bottom navigation bar has three icons: 'Record' (a blue circle), 'Sessions' (a list icon), and 'Settings' (a gear icon).

5:07    


[Edit](#)

Record Session

RINGS

[Add Senstream Ring](#)




SESSION DETAILS

Session Type Scoring: PPG 200Hz, EDA 5Hz, Temperature 5Hz 

Session Name

Notes

[Add a Senstream Ring to Continue](#)

Record Sessions Settings

Settings Page

