

	Labelling	Référence
	<p>Plan Labelling</p> <p>e-Celsius® Médical - Capsule</p>	P022LAB003

Cartouche d’approbation du document

	Nom et prénom	Fonction	Date	Signature
Rédacteur	BLOND Estelle	Directrice Technique		
Vérificateur	SEYCHELLES Jean-François	RQAR		
Approbateur	BLOND Estelle	Directrice Technique		

Suivi des modifications du document

Révision	Date de modification	Description des modifications
1	21-09-2020	Création

Page laissée libre volontairement.

1. Objet du document et périmètre d'application

This document describes the label plan for the blister of the capsule product e-Celsius® Medical.

2. Documents et acronymes utilisés

A. Documents applicables

P022SEE001_P022-M capsule Labelling requirements

B. Documents de référence

PRO18: Planifier, Concevoir et Développer le produit.

C. Liste des acronymes et définitions



1:1

Label Information :

Dimension : 40*40 mm

Support : laize PGL du blister (AMCOR Flexibles SPS ; référence 4217000224)

Impression : ruban transfert Thermique « Qualité » W
(Allen France ; référence FN2245W)

AAAA-MM: Expiration year AAAA and month MM

AA-XX: Batch number AA is the year and XX is an incremental number

YYYY: Notify body number

Data Matrix include:

(01)GTIN Code : 3665810000076

(37)Quantity of product (N=1)

(10)Batch number (AA-XX)

(17)Manufacturing date (YYYY-MM)

Code: 01366581000007637N10AA-XX17YYYY-MM

Position on the product :

This label is printed on the capsule blister.

