



Quick-start guide

Welcome




Get Started

- 1 Go to **greatcall.com/wearable/activation** to activate your Lively service.
- 2 Once activation is complete, download the Lively App to your smartphone.
- 3 Through the Lively App you will be guided on how to pair your Lively Wearable with your smartphone.


Wearing your Lively Wearable

Insert the Lively Wearable Jewel from behind the Wristband or Pendant accessories. Fall Detection only works when worn around the neck using the Pendant accessory.



Using your Lively Wearable

In an emergency, press the button located on the face of the Jewel and you will be connected to a 5Star Agent on your smartphone.



For more information on how to use your Lively Wearable, you can refer to the User Guide online at greatcall.com/wearable/userguide.v



greatcall.com/support

(866) 490-1076