# PRODUCT INFORMATION

SPECIFICATIONS:

Company Name: PillDrill, Inc.

Name and Model: PillDrill Hub, PD0002 Power Source: DC 5V 1A Adapter

AC Power Supply: 100-240V, 47-63Hz, 0.25A rms

Operating Temp: 0°C - 40°C

USAGE: Use the included power adapter and USB Cable. For further information, see the User Guide.

#### COMPLIANCE:

USA: FCC and CANADA: IC

FCC/IC WARNING

This device complies with part 15 of the FCC Rules and RSS-247 of Industry Canada. Operation is subject to the following two conditions:(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiveris connected.
- Consult the dealer or an experienced radio/TV technician for help.

CAUTION: Any changes or modifications to this device not explicitly approved by manufacture could void your authority to operate this equipment.

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes:

(1) l'appareil ne doit pas produire de brouillage, et

(2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.

### SAR Declaration:

This equipment complies with FCC/IC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance of 20 cm between the radiator and your body. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Cet équipement est conforme aux limites d'exposition aux rayonnements de la IC établies pour unenvironnement non contrôé. Cet équipement doit être installé et fonctionner à au moins 20cm de distance d'un radiateur ou de votre corps.

### FU: CF

PillDrill, Inc., hereby declares that the radio equipment type PD0002 is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address: www.pilldrill.com/compliance

Frequency Range:

Connectivity/Wi-Fi: 2.4GHz 802.11 b/g/n, 2412MHz ~2462MHz

Max power for Wi-Fi: 16 dBm

## Applicable CE Marking:





### RF EXPOSURE:

The Maximum Permissible Exposure (MPE) level has been calculated based on a distance of 20 cm between the device and the human body. To maintain compliance with RF exposure requirement, maintain a 20cm distance between the device and human body when using the product.

RECYCLING: Ensure that your PillDrill Hub is recycled or disposed of in accordance with your local laws and regulations.

**PILLDRILL USER GUIDE:** Prior to using the PillDrill Hub, please review the information in the User Guide located at www.pilldrill.com/userguide

WARRANTY: This Limited Warranty applies to physical goods purchased from PillDrill, Inc. (the "Product"). This Limited Warranty covers any defects in material or workmanship under normal use during the Warranty Period, During the Warranty Period, PillDrill, Inc. will repair or replace, at no charge, Products or parts of a Product that prove defective because of improper material or workmanship, under normal use and maintenance. PillDrill, Inc. will either repair the Product at no charge, using new or refurbished replacement parts. The Warranty Period for Products purchased from PillDrill, Inc. is 365 days from the date of purchase. A replacement Product or part assumes the remaining warranty of the original Product. This Limited Warranty does not cover any problem that is caused by conditions, malfunctions or damage not resulting from defects in material or workmanship. To obtain warranty service, you must first contact PillDrill, Inc. to determine the problem and the most appropriate solution for you.

**LIMITATION OF LIABILITY:** In no event shall PillDrill, Inc. be liable to you for any damages resulting from your purchase or use of the Product or related PillDrill, Inc. services such as the PillDrill app or website. In no event shall PillDrill, Inc. be liable for any injury, loss, claim or any indirect or consequential damages including but not limited to, lost profits, revenue, savings, data or for any other damages based either in contract law, tort law or statutory strict liability or otherwise. Where the exclusion or limitation of liability is not allowed by law, our liability shall be limited to the maximum extent permitted by law and in no event shall any liability of PillDrill, Inc. exceed the amount paid by a customer for the Product or related services.

