







 YYYY-MM-DD  
**QTY** 1

## AURA™ Intracranial Pressure Monitoring Kit

For use with AURA™ Monitor T-101.

<b>REF</b> SP101-A	<b>SN</b> XXXXXXXXXX	<b>LOT</b> XXXXXXXXXX
 YYYY-MM-DD	20 %  85 % 50 kPa  119 kPa	0 °C  50°C



**R<sub>x</sub>** ONLY



STERILE EO



UDI: (01)01234567890123(17)1234567890

### Package Contents:

- AURA™ ICP Sensor
- Drill Bit with Collar
- Hex Wrench
- Scalp Retractor
- Antenna Holster



FCC ID: 2AJW6-BTAIS-01

Patent Pending

Inside contents intended for single use only.

Distributed By: Branchpoint Technologies, Inc. 1 Technology Drive., Suite I-811  
Irvine, CA 92618 USA+1 949 829 1868

PN: EMT AD-1412-001 r.C

Figure 1 - Sterile Packaging Label for AURA(TM) Sensor





YYYY-MM-DD

QTY 1

# AURA™ Intracranial Pressure Monitoring Kit

For use with the AURA™ Monitor T-101.

<b>REF</b> SP101-A	<b>SN</b> XXXXXXXXXX	<b>LOT</b> XXXXXXXXXX
YYYY-MM-DD	85 % 119 kPa 50 °C 20 % 50 kPa 0 °C	



R<sub>X</sub> ONLY



STERILE EO



UDI: (01)01234567890123(17)1234567890

## Package Contents:

- AURA™ ICP Sensor
- Drill Bit with Collar
- Hex Wrench
- Scalp Retractor
- Antenna Holster
- Patient Registration Card
- Patient Labels (5)
- Instructions for Use



FCC ID: 2AJW6-BTAIS-01

Patent Pending

Before opening the package, run in-package check for pre-implantation status of AURA™ ICP Sensor. Follow steps in AURA™ Monitor Kit Instructions for Use to complete in-package check.

Distributed By: Branchpoint Technologies, Inc. 1 Technology Drive., Suite I-811  
 Irvine, CA 92618 USA+1 949 829 1868  
 PN: EMT AD-1413-001 r. C

Figure 3 - AURA(TM) Sensor Sales Packaging Label



Figure 4 - AURA(TM) Sensor Sales Packaging showing Label location