

Product: Wireless Probe Type Ultrasound Scanner  
Model:CP-probe 2 type CL  
Built-in Battery: d.c.3.85V 2800mAh  
Power Energy: 5VA Rated Input power: d.c.4V,1A  
FCCID:2AKFH-CPROBE

 Guangzhou Sonostar Technologies Co., Ltd.  
504#, C Building, #27 Yayingshi Road, Science  
Town Guangzhou, CN 510665  
Tel: +86 20 32382095 Fax: +86 20 62614030  
Manufacturer Website: www.sonostar.net

SN  Check  
probe body   Check  
probe body  
Serial number Date of manufacture

**Prescription Statement**"Caution:Federal law  
restricts This device to sale by or on the  
order of a physician"

#### Operation

1. If use iPhone/iPad, search "WirelessUSG" on APStore, if by Android copy software to phone
2. Press key on probe, power on (long press for power off, short press for Freeze), go to Wifi setting of phone/tablet, click wifi from probe and input its code number as password.