

VisionAR

QUICK GUIDE

ETL-PN0005.00 Rev.1.2

What's in the box?


LEVEL 1:

- Elastic Headband
- Clip
- USB cable
- Notes
- Cleaning Cloth and liquid
- Screwdriver


LEVEL 2:

- Glasses with lens
- Spare lens
- Case


1 • Attach the headband




2 • Plug in the USB C




3 • Put your VisionAR on




4 • Adjust the Headband to secure your VisionAR



5 • Connect the magnetic clip



6 • Connect the CU and VisionAR



For more information scan here!

