











Device label

MD REF	Nerivio™ device Model FGD000075	 www.theranica.com	   
Theranica Bio-Electronics LTD. 4 Ha-Omanut St., Netanya 4250438, Israel FCC ID: 2AOH8-NM Battery: 1 x Li-MnO2 3.0V, 1.2Ah Battery Max. Voltage: 3.3V --- Device Max. current: 40mA Frequency: 100-120Hz <small>LBL-NM-0010-2.0</small>		IP22 CEpartner4U BV Esdoornlaan 13, 3951 DB Maarn, The Netherlands	

Device label with UDI example

MD REF	Nerivio™ device Model FGD000075	 www.theranica.com	   
Theranica Bio-Electronics LTD. 4 Ha-Omanut St., Netanya 4250438, Israel FCC ID: 2AOH8-NM Battery: 1 x Li-MnO2 3.0V, 1.2Ah Battery Max. Voltage: 3.3V --- Device Max. current: 40mA Frequency: 100-120Hz <small>LBL-NM-0010-2.0</small>		IP22 CEpartner4U BV Esdoornlaan 13, 3951 DB Maarn, The Netherlands	
		UDI GTIN (01) 0 7290017 578019 (17) 210226 SN (21) NM2003-0000049 MAC address: F8:F0:05:A6:33:52 Use by: 2021-02-26	

