



Quick Start Guide

INVENSOR

- 1 Install the Invensor App and click on Create Account. Enter your information and click done, the Invensor App will automatically begin scanning for any available Invensor units. Make sure that Bluetooth is enabled on your device and the Invensor unit is within range.



- 2 Once connected, the Invensor App will prompt you to enter information associated with the Invensor unit you connected to.




- Where will this Invensor be placed?
- Add medicines and set temp range for each.
- Set Blood Glucose range.
- Add Caregiver

- 3 Once complete, this process will repeat for each Invensor unit you are attempting to connect to your device.

Your device is now set up to track and manage your medication, blood glucose, reminders, alerts, and peace of mind.

Manage your medication.
Worry less. Live Life.



FCC

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna. Increase the separation between the equipment and receiver. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. Consult the dealer or an experienced radio/TV technician for help.

ISED

User Manual Notice for License-Exempt Radio Apparatus (RSS-Gen, Issue 4, Nov. 2014, Section 8.4)

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions:

1. This device may not cause harmful interference;
2. This device must accept any interference received, including interference that may cause undesired operation of the device.


Cet appareil est conforme à Industrie Canada une licence standard RSS exonérés (s). Son fonctionnement est soumis aux deux conditions suivantes:

1. Cet appareil ne doit pas provoquer d'interférences;
2. Cet appareil doit accepter toute interférence reçue, y compris les interférences pouvant provoquer un fonctionnement indésirable de l'appareil.

INVENSOR

Download the app

Before you use your Invensor you need to install the app and sync your device via Bluetooth. Please refer to the steps in the Quick Start Guide for more info. You can download the free App on the App Store or Google Play. Simply search for "Invensor" or click the associated links at www.invensify.com.



INVENSOR

Quick Start Guide

Manage your medication.
Worry less. Live life.