



11.c




Start manual mode and see the stimulation's pace

12




View patient's **history**

13



View patient's **performance**

14



View **details** for each exercise

WARNING

Make sure the device battery is fully charged before each session.

MySPM requires internet connection and active bluetooth.

MySPM requires confirmation of location and storage permissions.

MySPM

Shoulder Pacemaker™ App

REF **SPM01A1**



NCS Lab srl

Via Pola Esterna 4/12

41012 Carpi (MO) - Italy

Ph: +39 059 669813

info@ncs-company.com



0476

www.shoulderpacemaker.com

info@shoulderpacemaker.com