








Genrui Auto Hematology Analyzer
Model: KT-6810 [REF] 21000013
Power: 100-240V~50/60Hz 200VA
CE
YYYYMMDD
Genrui Biotech Inc.
Address: 10, Shenzhou Building B, Nanhai Tangle Innovation
Industrial Park, 1183 Lianhai Road, Nanhai District, 510555, Shenzhen, China

Class 1
Laser product

Power Input: 100-240V~50/60Hz
Fuse: Plug F6-3AL250V

HIGH VOLTAGE

First Startup Notice
1. Open the left door, cut the line leading to the
jet assembly of the reagent container.
2. Open the right door, cut the line leading to the
sampling module and remove the test
3. Remove the cartridge. See on the injection
4. Close the left and right doors.
5. After the first startup, turn on the power.

DILUENT

WASTE

