


Description of E-Label

FCC ID : 2BBVT-DSI-100

IC: 30869-DSI100

This device is applied E-label.

Step 1. Run scan app settings.



Step 2. Select the General information tab.

The screenshot shows a 'Settings' window with a sidebar on the left and a main content area on the right. The sidebar has a 'General' tab selected, which contains 'Information' and 'Notation' sub-items. The main content area is titled 'Information' and contains three sections: 'Version', 'Certification', and 'About'. Red circles with numbers 1, 2, and 3 are placed next to these sections respectively. The 'Version' section shows 'ScanAppi V0.0 (Build 0.00.00.0000)'. The 'Certification' section shows a URL 'https://doflab.com/certification/'. The 'About' section shows 'ScanAppi™', a URL 'http://www.doflab.com', and a copyright notice 'Copyright 2023. DOFlab Co.,Ltd. All rights reserved.'.

Settings

General | System | Scanner | Back up

General

Information

Notation

1 Version

ScanAppi V0.0 (Build 0.00.00.0000)

2 Certification

<https://doflab.com/certification/>







3 About

ScanAppi™
<http://www.doflab.com>
Copyright 2023. DOFlab Co.,Ltd. All rights reserved.






Sign out Save Cancel

Step 3. Go to web site link number 2. (<https://doflab.com/certification/>)

Certification

Symbol	Description
	Manufacturer and address
	Date of manufacture and country of manufacture (Korea)
	European representative and address
	Serial Number
	Medical device standard code
	Medical Equipment

The [open recovered files disk drill](#) author also highlighted the significance of regularly backing up files to prevent such mishaps in the future. Is it possible [how to recover deleted files from usb online](#) to recover files from a formatted flash drive on Windows 8?

South Korea	U.S	Canada	Europe	Japan
 모델명 :DSI-100 제조사:DOF INC 제조년월:별도표기 인증자상호:DOF INC 명칭:구강스캐너 제조국:한국 유선전화:070-XXX-XXXX R-C-dof-DSI-100	 FCC ID : 2BBVT-DSI-100	PMN : FREEDOM i HVIN : DSI-100 IC : 30869-DSI100	 	 020-230340

There are no special access codes, accessories, or permissions beyond the normal security protection to unlock the screen, sign-in page, or overall product access. The E-Label will be easily accessible by the user and cannot be modifiable by the user
The E-Label will be clearly legible without the aid of magnification.

This device can only be operated via a dedicated software.

Regulatory information can be checked through the website linked to the dedicated software.

Regulatory information cannot be changed or deleted except the manufacturer.

BOX Label

Product name : FREEDOM i

Model name : DSI-100

Rating

Input : (100-240) V~, (50/60) Hz, 1.2 A

Output : 12 Vd.c., 3.4 A

Intended use

A device that takes a video image of the condition of the oral cavity using a camera and reproduces it in a three-dimensional image.

EC

REP

DOF Inc #B204, 601, 602, 77, Seongsuil-ro, Seongdong-gu, Seoul, Republic of Korea

+82 70-5057-3518

www.doflab.com

FCC

CE

FCC ID : 2BBVT-DSI-100

XXXX

R-C-dof-DSI-100

IC : 30869-DSI100

MD

Rx ONLY

NON STERILE

-5C

800hpa

45C

1100hpa

20%

80%

xxxx-xx-xx

BE	BG	CZ	DK
DE	EE	IE	EL
ES	FR	HR	IT
CY	LV	LT	LU
HU	MT	NL	AT
PL	PT	RO	SI
SK	FI	SE	NO
IS	LI	CH	TR

QR 코드라벨

Intraoral Scanner