



Global Product Certification  
EMC-EMF Safety Approvals

## EMC Technologies Pty. Ltd.

ABN 82 057 105 549

### Melbourne

176 Harrick Road  
Keilor Park, Vic 3042  
Tel: +61 3 9365 1000

### Sydney

Unit 3/87 Station Road  
Seven Hills, NSW 2147  
Tel: +61 2 9624 2777

Email: [emc-general@emctech.com.au](mailto:emc-general@emctech.com.au)

Web: [www.emctech.com.au](http://www.emctech.com.au)

<b>To</b>	Federal Communications Commission
<b>Attention:</b>	Certification Officer
<b>EMCT Ref:</b>	M2211017
<b>Date:</b>	September 24, 2024
<b>Subject:</b>	Tournicare Pty Ltd Request for Original Certification - FCC ID:2BD2B-ARMA-08

<b>Product Name:</b>	ARMA Blood Pressure Monitor
<b>Model:</b>	08
<b>FCC ID:</b>	2BD2B-ARMA-08
<b>Transmitter:</b>	STM32WB55RGV6 Bluetooth Low Energy (BLE v5.2)
<b>Operating band:</b>	2400 – 2483.5 MHz
<b>Antenna:</b>	AN3359 copper trace (STMicroelectronics) – 1.96 dBi
<b>Description: (supplied by customer)</b>	The ARMA device is a blood pressure monitor designed for home use only. Results are shown on the LCD display as well as sent via Bluetooth to a connected smart phone once the test is complete.

### Please Note:

<b>Permanent Confidentiality:</b>	Application includes request for permanent confidentiality – Operational Description, Schematics, Block Diagram
<b>Short-Term Confidentiality:</b>	Application includes request for short-term confidentiality (180 days)  Internal photos: CONFE2207-1563-3 Annex B.pdf User Manual: ARMA-IFU, Version 4.pdf
<b>FCC 15.19 statement – exclusion request:</b>	Application includes request to omit the FCC 15.19 (a)(3) statement from the device label.
<b>Agent:</b>	EMC Technologies is acting as agent for Tournicare Pty Ltd in this application.

Regards

Shabbir Ahmed  
CEO  
EMC Technologies Pty Ltd  
176 Harrick Road, Keilor Park, VIC 3042 Australia  
Ph: +61 3 9365 1000 [www.emctech.com.au](http://www.emctech.com.au)