

# ARMA LABEL



REF

ARMA-08

Automatic Blood Pressure Monitor



Tournicare Pty Ltd  
18 Elliot Avenue, Balwyn  
Victoria, Australia 3101



DD MMM YYYY

FRM-063, Version 3



MD

Type BF Applied Part  
Rating: IP22 DC 5V 500mA  
FCC ID: 2BD2B-ARMA-08

SN

UDI

PI 0086274-C

LABEL MATERIAL: THERMLfilm SELECT 22720  
MANUFACTUROR: Flexcon  
LABELLING METHOD: The labelling method is a self adhesive label. Using a thermal transfer print method.



LABEL POSITIONED ON THE INNER ASPECT OF THE LARGER LEG OF THE C-CLAMP.



**REF** ARMA-08

Automatic Blood Pressure Monitor



Tournicare Pty Ltd  
18 Elliot Avenue, Balwyn  
Victoria, Australia 3101



DD MMM YYYY

FRM-063, Version 3



**MD**

Type BF Applied Part  
Rating: IP22 DC 5V 500mA  
FCC ID: 2BD2B-ARMA-08

**SN**

**UDI**

PI 0086274-C