


Label inside the shoe (for all shoes)

magnes


- Stitched
- Original size 3.5 x 2.5 cm



Label on charging station


Label on charging station

- Lazer engraved
- Original size 12.9 x 2.6 cm



Label on power adapter


- Lazer engraved
- Original size 3.8 x 2.2 cm



Label on box - NUSHU X

magnes


- Sticker
- Original size 9.4 x 7 cm



Label on box - NUSHU FDA

magnes

- Sticker
- Original size 9.4 x 7 cm



Label on box - NUSHU CE MD

magnes

- Sticker
- Original size 9.4 x 7 cm

NUSHU

Shoe for gait data collection

MD



SN

 XX XXXX

REF

 010.001.010

CE

 0297

UDI



(01)xxxxxxxxxxxxxxxx(11)xxxxxx(21)xxxxxxxxxx

010.001.007.v1

 IP22  1500 mbar 70%
600 mbar 0%   Max 120 kg  0°C +30°C  3.6 V
Li-ion 750 mAh 

 Effectum Medical AG
Kirchgasse 11
CH-4600 Olten

EC

REP

MED-RAS GmbH
Eichenallee 8H
DE-21521 Wohltorf

magnes
www.magnes.ch