

Date: 2024-08-05

TUV SUD BABT UNLIMITED TCB  
Octagon House,  
Segensworth North,  
Fareham, Hampshire,  
PO15 5RL, United Kingdom

We, ShenZhen GoodlyMed Technology Co.,Ltd. (Address: 701, Building C, Area C, Datianyang Industrial Zone, Shiwei Community, Matian Street, Guangming District, 518107 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA )

Hereby declare:

Product name: Upper Arm Electronic Blood Pressure Monitor  
FCC ID: 2BERFA01SLB001


All models are identical in RF module, PCB layout, main schematic, critical components and main software. The differences are listed in the below table.

**Model list for below table:**

Model name	A01-SLB	A01-BS	A01-BL
Dimension	129x97x46mm	129x97x46mm	129x97x46mm
Power supply	The same		
Voice	With	With	Without
Backlight	with	Without	with
Bluetooth	With	With	With
Hardware	The same		
Schematic	The same		
PCB layout	The same		
Software	Sampling circuits, algorithms are the same. The difference for models A01-SLB, A01-BS and A01-BL are only in the voice and backlight functions.		

Yours sincerely,

Company Officer: Maria Ma



Telephone Number: 18002580887

Email: mary@goodlymed.com