

Sartorius Stedim Systems GmbH, 34302 Guxhagen

Anton Matrosov  
Anton.Matrosov1@sartorius.com

11. November 2024

Sartorius Biowelder S - ACCESS AND INSTRUCTIONS TO DISPLAY  
REGULATORY INFORMATION.

As demanded by the FCC according to chapter 4.1 of the electronic labeling  
guidance:

#### **4.1 User Steps to Access Information**

Users must be able to access the regulatory information without requiring  
special access codes or permissions, and in all cases the information must be  
accessible in no more than three steps from a product's main or home menu.  
An example access process is as follows.

#### **Starting from Main or Home Menu:**

Step 1 – accessing submenu (for example, "Settings")

Step 2 – access submenu (for example, "Legal and Regulatory")

Step 3 – access submenu (for example, "Regulatory", where applicable  
regulatory information is displayed)

Follows the photo documentation.

Best Regards  
Sartorius Stedim Systems GmbH

Anton Matrosov

Sartorius Stedim Systems GmbH  
Robert-Bosch-Straße 5-7  
34302 Guxhagen, Germany  
Telefon +49 5665 407 0  
Fax +49 5665 407 2200  
sartorius.com

Sitz der Gesellschaft:  
Guxhagen  
Registergericht:  
Amtsgericht Fritzlar  
HRB-Nr. 11105

Geschäftsführer:  
Ralf Slowik,  
Jens Karl

Vorsitzende des Aufsichtsrates:  
Michaela Pischke


Step 1: Open Menu by pushing on (1)

Blade count: 5/50


## Tube welding



Step 2: Open „Device Informations“ by pushing on (2).



Step 3: Scroll down to Section “Regulatory Information” and pushing on them (3).



Display of “Regulatory Information’s”:

