



# Device pack container label

Placed on the pouch that the device is packaged in.




# Device box label

Placed on the outer box that the device is packaged in.


**Perin Health**

Perin Health Patch



Perin Health Devices, LLC  
755 Research Parkway, Suite 530  
Oklahoma City, OK 73104  
[www.perinhealth.com](http://www.perinhealth.com)


SN

  
A1234567890A


REF


PHP80060


UDI





(01)01234567891011  
(17)210526  
(10)123456789


  
Single Patient Use Only


  
Single Patient Use Only








  
LATEX

  
Keep Dry

  
-20 °C +50 °C





  
eIFU indicator

Contains 1 x Perin Health Patch  
(Rev. A | 80060-2)

IP44

Made in the USA

This device complies with Part 15 of FCC rules. Operation is subject to two conditions: 1) this device may not cause harmful interference, and (2) this device must accept any interference that may cause undesired operation.  
**FCC ID: 2BKSW-PHP-0124**

120 mm

50 mm