

MEMTEC CORP

INSTRUCTIONS

- PREPARE THE PATIENT FIRST
- CLOSELY FOLLOW INSTRUCTION CARD IN THE KIT
- INSTALL FLASH CARD WITH LABEL FACING UP
- INSTALL A NEW OR RECHARGEABLE 1.5 VOLT AA BATTERY
- REVIEW LCD STATUS
- LOG DIARY TIME AND PLACE HOLTER PERFORMER IN POUCH

MEMTEC CORPORATION

68 STILES RD, SALEM, NH 03079

PH (603) 893-8080

OPEN

