



MINI MITTER

20300 Empire Avenue, Building B-3 • Bend, OR 97701 USA • Fax (541) 322-7277 • Phone (541) 322-7272
email: mm@minimitter.com • web page: <http://www.minimitter.com>

October 6, 2003

Attention FCC Application Reviewer:

To comply with 47 CFR 2.925, paragraph (a), it is required that the transmitting device be labeled with an FCC Identifier which as quoted from subparagraph (1) must be "...preceded by the term FCC ID in capital letters on a single line, and shall be of a type size large enough to be legible without the aid of magnification."

Our device, the ingestible core temperature capsule, is only 22mm x 8mm ϕ . Therefore, it is impractical to label this device with the FCC ID. Additionally, the capsule transmitter is intended to be ingested, and therefore, the FCC ID would not be visible when in use. Following the directions stated in paragraph (e), Mini Mitter presents an alternative to a permanently affixed nameplate. As an alternative, Mini Mitter proposes to place the required FCC ID on the envelope in which the capsule transmitter is packaged. The capsule remains in this envelope until ready for use, at which time it is activated and swallowed. Please see the drawing for the proposed capsule label using the link below. Also below are photographs of the capsule packaging.

With best regards,

Donna Barton
RF Design Engineer
Mini Mitter Co., Inc.
20300 Empire Ave. Bldg B-3
Bend, Oregon 97701
<http://www.minimitter.com>

Attachments.

Link to Drawing of Proposed Capsule Label showing FCC ID#


[Capsule Label.pdf](#)



MINI MITTER

20300 Empire Avenue, Building B-3 • Bend, OR 97701 USA • Fax (541) 322-7277 • Phone (541) 322-7272
email: mm@minimitter.com • web page: <http://www.minimitter.com>

Figure 2: Photograph of Proposed Capsule Envelope Label showing FCC ID#. The FCC label ID# is on the row third from the bottom. A photograph of the capsule sensor assembly and subassembly is placed next to the envelope for size comparison.





MINI MITTER

20300 Empire Avenue, Building B-3 • Bend, OR 97701 USA • Fax (541) 322-7277 • Phone (541) 322-7272
email: mm@minimitter.com • web page: <http://www.minimitter.com>