



Rayence Co.,Ltd.

14, Samsung 1-ro 1-gil, Hwaseong-si, Gyeonggi-do, Korea
Tel. +82-31-8015-6459 / Fax +82-31-8015-6394

Date: November 8, 2016

To: SIEMIC, INC.
775 Montague Expressway,
Milpitas, CA 95035
USA

Dear Sir/Madam,

Re: Letter of Similarity

We, Rayence Co.,Ltd., declare the differences for between original model and variant models.

-Product: Wireless digital X-ray flat panel detector

-Model name: 1717WCC

-Variant Model name: 1717WGC

There's different Sintillator type.

1717 WCC (Basic M/N) & 1717 WGC

CsI(Cesium Iodide/1717 WCC) or Gd₂O₂S:Tb(Gadolinium Oxysulfide/1717 WGC)

Sincerely Yours,

A handwritten signature in black ink, appearing to be "Ha" or "Namyoung Ha", written in a cursive style.

NAME : Namyoung, Ha

TITLE : Assistant Manager

E-MAIL: namyoung.ha@rayence.com

FAX: + 82-31-8015-6394

TEL: + 82-31-8015-6459

Rayence Co., Ltd.