

## 1 PURPOSE

This document is part of the Cochlear radio certification activities and shows external photos of the SM-2P.

### Sales Item numbers

Cochlear Wireless Mini Microphone 2+ GB	P770845
Cochlear Wireless Mini Microphone 2+ EU	P770846
Cochlear Wireless Mini Microphone 2+ AUS	P770847
Cochlear Wireless Mini Microphone 2+ US	P770848



Back



Bottom view (Left: line-in - right: FM/EURO Pin connector)  
LED ON .



Top view – Charging  
LED ON



Right side view. (Micro USB, Volume down, Mute button/LED, Volume up)



Left side view.

Pairing button  
(white)

Mode button  
(grey)




Mode LEDs (top – down):

- Microphone
- Telecoil
- Line-in
- FM

Back side, bottom end.



Mic LED ON



Telecoil LED ON



Line-in LED ON



FM LED ON



Back side – Clip design