



12151 Research Parkway, Suite 200  
Orlando, FL 32826 USA

P +1-800-331-3220  
P +1-407-275-3220  
F +1-407-249-2022


[www.invivocorp.com](http://www.invivocorp.com)

## Location of FCC Label


November 22, 2010

To Whom It May Concern:

The label bearing the FCC and IC identifiers will be located on the Essential MRI Patient Monitor as indicated on the following photograph, on the rear of the device beneath the battery during normal use:



Although the photograph shows a label that does not bear the FCC and IC identifiers, the production devices will bear the FCC and IC identifiers on this label as indicated in the following drawing:



## Location of FCC Label

The location of the FCC and IC identification numbers is consistent with the regulations which state that the identifiers be readily visible and that visibility at all times during normal use is preferable but not a prerequisite for grant of equipment authorization.

Invivo currently has several other wireless devices with market clearance with this identical structure: FCC and IC identifiers on a label separate from the model label, but located near or adjacent to the model label. Due to space constraints, there is no way to put all of the information on the same label because the device is not large enough to bear a label of such size. Therefore, the FCC and IC identifiers must be located on a label beneath the user-replaceable battery. Further, the following design specifications ensure that the FCC / IC identifiers will be visible by the operator often during device use:

- The minimum battery life is eight hours and the battery must be removed to be charged so the operator is guaranteed to see the FCC and IC identifiers at least two times per day, at the end of the day when removing for charging and the beginning of the day when inserting in the device.
- The device is handheld.
- The FCC and IC identifiers are visible at the time of purchase because the device is shipped to the user without the battery installed. (Because the battery is lithium-ion, Invivo is legally required to ship the device without the battery charged or installed, and the battery is required to be located in a box separated from the device.)
- The FCC and IC identifiers will be located on the device label and in the user's manual.
- The FCC and IC identifiers are located on a label which is not on a removable part. The FCC and IC identifiers are visible when the user-replaceable battery is removed and the battery is removable without the use of a tool.

*Krystal Mitchell*

Krystal Mitchell  
Quality & Regulatory Engineer  
Invivo Corporation  
a division of Philips Medical Systems