

DENTSPY Cavitron® Plus

Ultrasonic Scaling Unit

Type GEN.-136



Consult Instructions
for Use

100-240V 50/60Hz 1.0A 80VA
U/S POWER 30 WATTS
U/S FREQ. 30KHz

This device complies with Part 15
of the FCC Rules. Operation is
subject to the following two
conditions:

- 1) this device may not cause
harmful interference, and
- 2) this device must accept any
interference received, including
interference that may cause
undesired operation.

FCC ID: TF3-DPD81842
IC: 46818-DPD81842
Model: G136



MEDICAL EQUIPMENT
UL 60601-1

CAN/CSA-C22.2 No. 601.1
ANSI/AAMI ES60601-1 (2005, 3rd ed.)
CAN/CSA-C22.2 No. 60601-1 (2008)
13VA



AC Power



Type B Equipment



Wireless Device



Footswitch IPX-1



Manufactured by:
DENTSPY Professional
1301 Smile Way
York, PA 17404



DENTSPY DeTrey GmbH
De-Trey-Str.1
78467 Konstanz
Germany



WEEE

P/N 81857 Rev. 6 (10/13)

DENTSPLY Cavitron® JET Plus

Ultrasonic Scaling and Air Polishing Unit

Type GEN.-137



Consult Instructions
for Use

100-240V 50/60Hz 1.0A 80VA
U/S POWER 30 WATTS
U/S FREQ. 30KHz

This device complies with Part 15
of the FCC Rules. Operation is
subject to the following two
conditions:

- 1) this device may not cause
harmful interference, and
- 2) this device must accept any
interference received, including
interference that may cause
undesired operation.

FCC ID: TF3-DPD81842

IC: 46818-DPD81842

Model: G137



MEDICAL EQUIPMENT
UL 60601-1

CAN/CSA-C22.2 No. 601.1

ANSI/AAMI ES60601-1 (2005, 3rd ed.)

CAN/CSA-C22.2 No. 60601-1 (2008)

13VA



AC Power



Type B Equipment



Wireless Device



Footswitch IPX-1



Manufactured by:
DENTSPLY Professional
DENTSPLY International
1301 Smile Way
York, PA 17404



DENTSPLY DeTrey GmbH
De-Trey-Str.1
78467 Konstanz
Germany



WEEE

P/N 81883 Rev. 6 (10/13)

DENTSPLY Cavitron® Prophy-JET®

Air Polishing Unit

Type GEN.-138



Consult Instructions
for Use

100-240V 50/60Hz 1.0A 80VA

This device complies with Part 15
of the FCC Rules. Operation is
subject to the following two
conditions:

- 1) this device may not cause
harmful interference, and
- 2) this device must accept any
interference received, including
interference that may cause
undesired operation.

FCC ID: TF3-DPD81842

IC: 46818-DPD81842

Model: G138



MEDICAL EQUIPMENT

UL 60601-1

CAN/CSA-C22.2 No. 601.1

ANSI/AAMI ES60601-1 (2005, 3rd ed.)

CAN/CSA-C22.2 No. 60601-1 (2008)
13VA



Wireless Device



AC Power



Type B Equipment



Footswitch IPX-1



Manufactured by:
DENTSPLY Professional
DENTSPLY International
1301 Smile Way
York, PA 17404



DENTSPLY DeTrey GmbH
De-Trey-Str.1
78467 Konstanz
Germany



WEEE

P/N 81958 Rev. 4 (10/13)



DENTSPLY Cavitron® JET Plus
Ultrasonic Scaling and Air Polishing Unit
Type GEN.-137

12



Caution, Consult
Accompanying Documents

100-240V 50/60Hz 1.0A 80VA
U/S POWER 30 WATTS
U/S FREQ. 30KHz

This device complies with Part 25
of the FCC Rules. Operation is
subject to the following two
conditions:

1) this device may not cause
harmful interference, and
2) this device must accept any
interference received, including
interference that may cause
undesired operation.
FCC ID: TTS-DPO73227323
IC: 46818-73227323

FC



MEDICAL EQUIPMENT
(UL60601-1 CAN/CSA C22.2 NO. 60601
130A)



Wireless
Device



AC Power



Type B Equipment



Footswitch SPN-1

Manufactured by:
DENTSPLY Professional
DENTSPLY International
York, PA 17404
CE 0086

EC REP

DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
Germany

P/N 81803 Rev. 1 (10/11)

Made in U.S.A.



WEST