

D-Tag

Models

T33-110 – TeleMonB Telemetry Transmitter

Installation Instructions



RadarFind™ Corporation
877-RadarFind (723-2734)
www.radarfind.com



Warnings:

It is important to follow all instructions shipped with this product.

UNPACKING: After unpacking the unit, examine it for damage that may have occurred in transit. If the equipment has been damaged, do not attempt to install or operate it. File a claim immediately with the carrier stating the extent of the damage, or contact RadarFind customer service for repair or replacement.

Read and understand all instructions before installing or operating this equipment.

The selection of the mounting location for the device is to be accomplished under the direction of the director of Clinical Engineering and conform to all company policies, and must be in accordance with specifications in the sub-chapter for each device type.

CAUTION: Risk of malfunction. Install these devices only in areas with an ambient temperature between 0° and 50°C.

[15.19] This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

[15.21] Information to user: Changes or modifications not expressly approved by the RadarFind Corporation could void the user's authority to operate the equipment

[15.105] NOTE: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

[15.249] This equipment operates under Section 15.249 of the Title 47 FCC Rules.

This product protected by one or more of the following U.S. Patents: 7372365, 7525432 & 7528716. Other patents pending.

Contents of standard package:



Telemetry tag with stainless steel mounting bracket attached.

Accessories



TelemonB Installation Tool

General Guidelines:

Using stainless steel mounting bracket:

Step # 1 - Align the tag as shown in the picture and insert installation tool with machined side toward tag. This will hold the tag bracket aligned with the leads collar.



Use a rubber or plastic mallet to seat the bracket firmly onto the leads collar.

Step # 2 – (Assuming that you have already entered your assets - See next page for *ADDING* assets in the RadarFind database.) Select asset by entering the serial number or asset inventory number. Enter or scan tag ID into the Radarfind software, select which tag type, status or non-status, and hit <Tag!>

hardware > assets

Asset Type: Asset SubCategory: ☐ Tagged Only ☒ Un-Tagged Only

Delete selected **Add an Asset** **ASSETS**

<input type="checkbox"/>	Name	Description	Type	Inventory(CE)Id	Manufacturer	Serial	Status	Tag	Where	Actions
<input type="checkbox"/>	BED		General Bed	FNB2701		1130-62A30-A	<input type="radio"/> With Status <input type="radio"/> No Status	<div style="border: 1px solid black; padding: 2px;">Tag!</div>		EDIT

Handwritten red text: TAG! with an arrow pointing to the Tag! button.

Once a tag is assigned to an asset any searching or display can be done by looking for the asset.

If you need to remove a tag from an asset please un-assign the asset in the following screen by selecting the correct asset and clicking DETACH.



▼User Interface | ▼Reporting | ▼Facility | ▼Hardware | ▼System |

hardware > assets

Asset Type: Asset SubCategory: ☐ Tagged Only ☐ Un-Tagged Only

Delete selected **Add an Asset** **ASSETS**

<input type="checkbox"/>	Name	Description	Type	Inventory(CE)Id	Manufacturer	Serial	Status	Type	Tag	Where	Actions
<input type="checkbox"/>	4100 PCA	INFUSION PUMP, PAT.-CTRL. ANALGESIC-16-924A	PCA Pump	02739	ABBOTT LABORATORIES DIAG DIV	21529	Available	AS1	00:00:00:03:10:00:02:e6	Radio	EDIT Detach

Run-Time Asset Additions or Edits

To edit a current asset locate the asset in the asset table and hit the edit icon.

To add a new asset to the database select Add an Asset from the menu bar.



Selecting either Edit or Add will bring you to the following asset edit screen where you can enter or edit details concerning the asset. The important lines are:

- 1 – Display Name - The Asset Name
- 2 – Inventory ID – Your property tag ID number - this number must be unique to all assets.
- 3 – Type – Select the appropriate asset sub-category from the drop down list.

See Appendix A for more information on Asset Categories and Sub-Categories

If you wish you can also scan a tag ID into the Tag field here as well as the list screen above.

Please remember to select that this is a NON-STATUS tag.

Edit Asset - 4100 PCA

Display Name:

Description:

Serial:

Inventory ID:

Class:

Manufacturer Id:

Type:

Tag:

Clinical Zone:

The Clinical zone field indicates the department which “Owns” or is “Home” to the asset. At this time it is not used elsewhere in the RadarFind system.

Technical Specifications:**Models T33-120 Telemetry Transmitter Tag**

DC power:	Internal sealed lithium battery
Power consumption:	Average 15uA
Communications	902-928 MHz FSK
Transmit Power	+10dBm maximum - *nominal* TX field level is 70dBμV/m @ 1m
Antenna:	Unity gain isotropic
Communications:	FCC Subpart 15.249
Operating:	5 to 35° C, 5% to 90% RH, non condensing
Packaging:	ABS Plastic (UL Flammability Rating V0)
Accessories:	Adhesive pad, holding base, plastic mounting straps, and steel mounting straps
Dimensions:	Optional holding base (2.6"L x 1.5"W x 0.375"D) Tag T33-110 (2.75" L x 1.75" W x .450" D)

*FCC Label***RadarFind Corporation™**

Appendix A – Asset Type Categories

Organization of Assets into a logical structure of Categories is important to allow users to easily locate Assets they are searching for and get the most benefit from the RadarFind system.

Each RadarFind system contains a default Asset Type Category Tree database installed as an example of how Assets may be organized to make the system easy to use.

Note: Assets imported from existing Asset management systems maintained by Biomed or Materials will be imported into Categories that match the original system. In order for general users to be able to locate items some alteration in the Asset naming conventions used by BioMed may be required.

The Asset tree may be modified to suit specific requirements, within some general guidelines to assure proper functioning of the Asset search and utilization functions of the system.

The Asset Tree View Screen

Assets: Asset Categories | Asset Categories/Assets Tree | Asset Pools

hardware > Asset Categories and Assets

Asset Categories [ROOT]

- Anesthesia Related or Respiratory
- Beds
- Blood Warmers [1]
- Cardiology Devices
- Computers/AV
- EKG/Defibrillation Devices [51]
- Feeding Pumps
- Heating/Cooling Devices [23]
- Infusion Pumps [490]
 - Blood Pump
 - Breast Pump
 - Multi-Channel
 - Peristaltic Pump
 - Pumps, Pressure Pad
 - Single IV Pump [1021]
 - Syringe Pump
 - Thermobac Drainage Pump
 - Triple IV Pump [82]
 - Volumetric Infusion Pump
 - Miscellaneous Equipment [3]
- M
- M
- OR Related
- PCA
- Portable Radiology
- PT Related
- radarfind tag [8]
- Rental Equipment [81]
- Respiratory Devices [53]
- Scales/Lifts [25]
- SCD/CPM Devices
- Telecommunications [2]

Assets ☐ Tagged Only ☐ Un-Tagged Only

Delete selected Add on Asset

Name	Type	Inventory(CE)Id	Serial	Status	Type	Tag	Around	Actions
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08015	14124907CC	In Use	AS1	00:02:10:63	1PCH16	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08065	14125183CC	In Use	AS1	00:02:20:06	5MT1F	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08080	14124908CC	Dirty	AS1	00:02:07:06	3P46	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08142	14125184CC	Dirty	AS1	00:02:03:85	6N54	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08456	14124909CC	In Use	AS1	00:02:1C:91	2M31	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08415	14125185CC	Ready	AS1	00:02:17:03	4P116	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08172	14124910CC	Ready	AS1	00:02:1C:86	3PC01	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08060	14125186CC	Ready	AS1	00:02:1C:02	CNPT01	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08042	14124911CC	Dirty	AS1	00:02:1C:90	5N61	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08128	14125189CC	Ready	AS1	00:02:21:27	1PCH17	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08128	14125189CC	In Use	AS1	00:02:20:F6	2PC04	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08326	14125191CC	In Use	AS1	00:02:10:34	1P108	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08427	14124914CC	In Use	AS1	00:02:10:58	5N56	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08167	14125191CC	In Use	AS1	00:02:1C:CA	5N65	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08472	14124917CC	In Use	AS1	00:02:20:6C	4M14	Detach ChangeTag
<input type="checkbox"/> Single Infusion Pump	Single IV Pump	D08218	14125193CC	Non Status	AS1	00:02:1D:90	8MEP12	Detach ChangeTag

Found 487 results. Showing items 1 to 20
1 2 3 4 5 6 7 8 9 10 next set >>

Items / page: 20 change

RadarFind Corporation 2008

Asset Category Tree Basic Guidelines:

- 1 – The tree is two levels deep, a Major Category and a Sub-Category.
- 2 – An Asset Sub-Category must be assigned to a Major Category.
- 3 - Individual Assets **may only be assigned to an Asset Sub-Category**, not to a Major Category
- 4 - Asset Sub-Categories may be moved to other Major Categories by simply selecting another "Parent" in the dropdown box at the bottom of the Edit screen.
- 5 – Assets may be tagged directly on the Asset Category screen if desired.

Important Controls in the Asset Tree View screen

The screenshot displays the 'Assets, Asset Categories, Asset Categories/Assets Tree, Asset Pools' interface. The left pane shows a hierarchical tree of asset categories, with 'Single IV Pump [407]' selected under 'Infusion Pumps [490]'. The right pane shows the 'Edit Asset Category - Single IV Pump' form. Annotations with arrows point to specific controls:

- Edit Icon:** Points to the edit icon (pencil) next to the 'Single IV Pump' category in the tree.
- Add an Asset Icon:** Points to the 'Add' icon (plus sign) next to the 'Single IV Pump' category in the tree.
- Add Sub-Category Icon:** Points to the 'Add' icon (plus sign) next to the 'Infusion Pumps' category in the tree.
- Turn on Category in Views:** Points to the 'Active' checkbox in the 'Edit Asset Category' form, which is currently checked.
- Assign Sub to Major Here:** Points to the 'Parent' dropdown menu in the 'Edit Asset Category' form, which is currently set to 'Infusion Pumps'.

The 'Edit Asset Category - Single IV Pump' form includes the following fields:

- Name: Single IV Pump
- Description: Single Channel Infusion Pump
- Fast Access Code: [Empty]
- Max accepted utilization (%): 100.0
- Active: Yes (checked)
- MapView Viewable: Yes (checked)
- QuickView Viewable: Yes (checked)
- QuickViewCS Viewable: No (checked)
- TransportView Viewable: No (checked)
- Synonyms: [Empty]
- Separate by 'T': [Empty]
- Parent: Infusion Pumps

Buttons at the bottom of the form include 'Save asset type', 'Save And Merge', and 'Cancel'.

RadarFind Corporation 2009