

Model Differences


The Rapid Test Reader can distinguish the differences of colors from various biochemical reactions and analyze the reactants to provide qualitative, semi-quantitative and quantitative clinical diagnosis.

- Main model: CHR-631W
- The series model is defined as: CHR-63xy, x can be 0~9 for different shell colors, y can be W for enabling WiFi and bluetooth functions, or not shown for disabling wireless functions.

Model No. Difference item	Main Model	Series Model								
	CHR-631W	CHR-63xy (x can be 0~9 for different shell colors, y can be W for enabling WiFi and bluetooth functions, or not shown for disabling wireless functions.)								
		y = “W” (Example: CHR-633W)	y = empty (Example: CHR-631)							
PCB layout and the circuit diagram	o	o								
Components	o	o								
Material	o	o								
Function	o	o	Disabling the wireless function							
Shape & Color	o	Appearance difference in color.								
Other	N/A	N/A								
Note: (1) “o” means the item is same with Main model.										
(2) Color in series										
x	0	1	2	3	4	5	6	7	8	9
Upper shell	White	White	White	White	White	Black	Black	Black	Black	Black
Bottom shell	Golden	Silver	Blue	Green	Pink	Golden	Silver	Blue	Green	Pink

Tainan, Taiwan 2024/05/15

Issue place and date



Mei-Ying Sung
General Manager

Stamp

