

System description:

The HearIt All system consists of three main units:

Transmitter

Transmits sound from the television and/or conversations to the receiver via a built-in 2.4GHz wireless digital audio system (WDA)

Receiver with neck loop:

Receives sounds from the transmitter via both WDA and Bluetooth. It works with headphones or hearing aids in T-position (telecoil).

Charging stand:

Charging Stand is for receiver and transmitter charging. It transmits the signal from a connected landline to the receiver via Bluetooth.

HearIt All has four main functions:

- Amplified sound from television/radio
- Amplified conversation
- Amplified telephone conversation
- Amplified mobile telephone conversation