



For Veterinary Use Only!

Refer all service to the factory. No user
servicable parts. Use only specified power
supply/charger Use only 10.8V battery

(LI2025X-72C)



FCC ID: XMOIBEX
IC: 8512A-IBEX



Contains TX FCC ID: ED98LMX9838
Contains TX IC: 1520A-LMX9838

For Veterinary Use Only!

Refer all service to the factory. No user
servicable parts. Use only specified power
supply/charger Use only 10.8V battery

(LI2025X-72C)



FCC ID: XMOIBEX
IC: 8512A-IBEX




Contains TX FCC ID: ED98LMX9838
Contains TX IC: 1520A-LMX9838

Door is Not Removable

EI Medical Ibex Pro regulatory label artwork

PART NUMBER: **0000-00000_EIUS_IBEX FCC_LABEL_ART-R01**





Notes:

- Material:** Self adhesive silver mylar (Total thickness with adhesive no greater than 0.010")
3M 467 adhesive
- Colors:** All graphics and text to be PMS Black, screened % PMS White, or screened % PMS Black
- Font:** All text in outline form

EI Medical Ibex Pro regulatory label artwork

PART NUMBER: **0000-00000_EIUS_IBEX PROREGULATORY_LABEL_ART-R04**

Referencing Part Number: 0000-00000_EIUS_REGULATORY_LABEL-R02



Notes:

- Material:** *Self adhesive (no thicker than .015") white polycarbonate
3M 467 adhesive*
- Colors:** *Top color to be PMS Cool gray 9C
All graphics and text to be PMS Black, screened % PMS White, or screened % PMS Black*
- Font:** *All text in outline form*



This device complies with Part 15 of the FCC Rules subject to the following conditions:

- 1) This device must not cause interference, and;*
- 2) This device must accept all interference, including interference that may cause undesirable operation.*

WARNING:

Modification of this device without consent of the responsible party may void the users right to operate the device.