

This insert contains important information about your Intel-GE Care Innovations™ QuietCare system.

Pager Alert Updates

An updated Pager Alert Quick Reference (DHF-3736) is included in this kit. The information in the quick reference replaces the information shown in similar tables in the following guides:

- Intel-GE Care Innovations™ QuietCare® User Manual (DHF-3734)
- Intel-GE Care Innovations™ QuietCare® Training Guide for Resident Caregivers (DHF-3735)
- Intel-GE Care Innovations™ QuietCare® Resident Caregiver Training DVD (DHF-3754)
- Intel-GE Care Innovations™ QuietCare® Hardware Maintenance Guide (DHF-3756)

Updated versions of these guides will be sent to you when available.

FCC Part 15 User Information Update

Do not attempt to disassemble the Care Innovations™ QuietCare system or modify the devices in any way. Changes or modifications not expressly approved by Care Innovations could void the user's authority to operate the equipment.

This device complies with part 15 of the FCC Rules. Operation of this device is subject to the following two conditions:

- (1) This device must not cause harmful interference.
- (2) This device must accept any interference received including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.



WARNING: THE QUIETCARE SYSTEM IS NOT AN EMERGENCY RESPONSE OR ALARM SYSTEM AND IS NOT INTENDED FOR PROVIDING MEDICAL CARE.

Copyright © 2013 Intel-GE Care Innovations LLC. All rights reserved. Care Innovations, the Care Innovations logo, and the Caring Icon logo are trademarks of Intel-GE Care Innovations LLC. QuietCare is a registered trademark of Intel-GE Care Innovations LLC. Intel and the Intel corporate logo are trademarks of Intel Corporation in the United States and/or other countries, used under license. GE and the GE Monogram are trademarks of General Electric Company in the United States and/or other countries, used under license.

Manufacturer: Intel-GE Care Innovations LLC, 3721 Douglas Boulevard, Suite 100, Roseville, California 95661