

# ViroVac

by Buffalo Filter



**CLASS I EQUIPMENT.**  
**TYPE CF APPLIED PART.**

IPX1

**PROTECT AGAINST DRIPPING WATER.**



**ALTERNATING CURRENT.**



## DANGER

Risk of explosion if used in the presence of flammable anesthetics.



## CAUTION

Electrical shock hazard. Do not remove cover. Refer servicing to qualified service personnel.



## ATTENTION

**BEFORE CONNECTING READ INSTRUCTIONS**

### CLEANING INSTRUCTIONS:

**UNPLUG UNIT** - Wipe with damp cloth containing mild disinfectant solution, or soapy water. Wipe dry with a clean cloth. **DO NOT STEAM STERILIZE.**

**This device complies with Part 15 of the FCC Rules.** Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**Part:** VV120

**Serial:**

Mfg.

Date



Remote  
Activator



9D93

Medical Equipment with respect to electrical shock, fire and mechanical hazards only. In accordance with UL 60601-1 and CAN/CSA C22.2 NO. 601.1.



**Manufactured by:**  
**Buffalo Filter**  
595 Commerce Drive  
Buffalo, NY 14228  
1-800-343-2324 (USA)  
716-835-7000 (Int)  
[www.buffalofilter.com](http://www.buffalofilter.com)

100/120V~  
10 Amps  
50/60 Hz

## CAUTION

This contains a RFID wireless tag FCC ID: Y92-BF0002 and IC 9532A-BF0002.



## CAUTION

Grounding reliability can only be achieved when the equipment is connected to a receptacle marked "Hospital Grade" or "Hospital Only".

