



Customer:	LUNDERG
Date:	2021-03-25
Model:	WPAS-10-LD
Type:	Silk Screen
Art Ver.:	Art_WPAS-10-LD_SC (20210325)
Size:	10" x 30"
Color:	Pantone 309



CONTAINS TRANSMITTER MODULE.
OPEN POCKET TO ACCESS THE 3-WAY SWITCH AND SELECT THE PAD SETTINGS.

←

LED TX-only OFF Both-sides

→

LED TX-only OFF Both-sides

●

●

TX-only for king or queen bed

Both-sides for twin size beds

Do not turn the Transmitter OFF unless you know you won't be using it for an extended period of time.

If the system isn't working correctly, first ensure the monitor is ON and that the TX-only/OFF/Both-sides switch in this transmitter is NOT in the OFF position. There is a small red light next to the switch.

If this red light doesn't turn ON when you apply and release pressure on the pad, then the pad needs to be replaced. Remember, this sensor pad has a lifespan of 18 months, and batteries can't be replaced.

If you DO see a red light, the bed pad is working correctly. If experiencing other problems, clear the monitor's memory and pair the bed pad and the monitor again.

If you have any questions, contact us at questions@lunderg.com