



IP22



Product Name: Smart Blood Pressure Monitor

Model No.: ePA-08B4 DC 6V , 4 x AA batteries

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

SN: B20110302090100001



2021-03 FCC ID: 2AAEEPA08B4



Shenzhen Better Health Measurement and analysis Technology Co. LTD.
702/704, Bock C, Tsinghua Univ Science park, No. 13 Langhan Rd,
Nanshan District, Shenzhen, 518057, China



产品底部贴纸

尺寸: 49.5x20mm, 防水(膜在最上)

材质: 白底黑字, 无边框

备注: 产品批号B20110302090100001跳至B20110302090101000(印刷在膜下)

NO.	PROJECT TITLE 测试内容	DATE	APPROVAL 审核
01	按照GB18308-2013要求设计制造	20210307	
	按照GB18308-2013要求生产	20210308	

深圳市信步健康测量分析技术有限公司 SHENZHEN BETTER HEALTH MEASUREMENT AND ANALYSIS TECHNOLOGY CO., LTD.			
FILE NO.: 文件编号		REV: 版本号: 01.2	
TITLE: 标题: 吴山舞墨图		DESC: 描述:	
AUTHOR: 作者: 吴山舞墨			
DATE: 日期:			
MEASUREMENT: 测量尺寸: mm		DATE: 20210307 日期:	
MEASUREMENT BY: 测量者: mm		DATE: 20210307 日期:	
CHECKED BY: 审核者: mm		DATE: 20210307 日期:	
APPROVED BY: 批准者: mm		DATE: 20210307 日期:	
		TOLERANCE 公差	

