


Model No. : TechClip




View of Product-21 L ear




View of Product-22 L ear




*The test report is effective only with both signature and specialized stamp and The result(s) shown in this report refer only to the sample(s) tested. Without written approval of BlueAsia, this report can't be reproduced except in full. The results described in this report do not represent the quality or characteristics of the sampled batch, nor do they represent any similar or identical products that are not explicitly stated.*




View of Product-23 L ear




View of Product-24 L ear




View of Product-25 L ear



View of Product-26 L ear



View of Product-27 L ear




View of Product-28 L ear




View of Product-29 R ear




View of Product-30 R ear




View of Product-31 R ear




View of Product-32 R ear




View of Product-33 R ear




View of Product-34 R ear



View of Product-35 R ear



View of Product-36 R ear



View of Product-37 R ear