



ivWatch, LLC
1100 Exploration Way, Suite 209
Hampton, VA 23666

1.855.IVWATCH (489.2824)
Fax: 757.645.4760
www.ivwatch.com

Compatible Antennas

The only compatible antenna for the 13.56MHz transmitter in the ivWatch Model 400 Patient Monitor is a custom designed flexible antenna, ivWatch part number CE-1001000 (drawing attached).

Unique Antenna Connector

The RFID transmitter in the ivWatch Model 400 includes an antenna that is a custom-designed flexible circuit which connects to the circuit board by an FFC connector. The antenna is mounted on an internal surface in the Model 400's housing in a specific position, requires professional installation by an ivWatch technician.

Antenna Type

Inductive Loop Antenna

Antenna Size

0.795 x 0.400 in.