

70.0


55.0

Product name: Body Composition Analyzer
Product model: CP10N Product size: 471*449*865mm
Measurement Range: 2-200kg Power Supply: 3.7V 1200mAh
Service Life: 5 years lithium battery, DC 5V/1A
Production Date and Production lot number: See the label.
Note: See the instructions for other details
Registrant/Manufacturer: Shenzhen Yolanda Technology Co., LTD
Company Address: 201, Jinfulai Complex Building, No. 49-1, Dabao Road, Dalang Community, Xin'an Street, Bao'an District, Shenzhen
Manufacturing Address: Factory C, No. 18 Guangyao Industrial Area, Zhulongtian Road, Shuitian Community, Shiyao Street, Bao'an District, Shenzhen, China
After-sales service hotline: 400-998-2913

Do not pull the wire

Warnings:

- Please follow the user manual for proper usage.
- Do not attempt to disassemble the equipment on your own, as doing so may lead to electrical shock accidents.
- Do not use the equipment in conjunction with other medical electronic devices, such as pacemakers.
- Keep the equipment away from liquids.
- Avoid applying excessive force when handling the hand electrodes.



FC
FCC ID: 2ANDX-CP10N

This device complies with part 15 of the FCC Rules. Operation is subject to the condition that this device does not cause harmful interference (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

4-R1

