

KYMATI GmbH · Am Hochacker 5 · 85630 Grasbrunn · Germany

MARTIN GLÄNZER

to Whom it may concern

+49 89 5157599-60
+49 151 44025870
Martin.glaenger@kymati.com
www.kymati.com

page 1 of 2

Grasbrunn, May 26, 2023

SUBJECT: DECLARATION OF SIMILARITY KY-LOC 1D.02.01

Dear Ladies and Gentlemen,
to Whom it may concern:

Please be advised that the Kymati sensors:

Product Marketing Name	Hardware Version Identification Number
KY-LOC 1D.01.01	1.1.1.1.2
KY-LOC 1D.02.01	1.1.1.1.2
KY-RAY 1D.01.02	1.1.1.1.2
KY-RAY 1D.02.02	1.1.1.1.2

are electrical identical. The difference between the sensors is the Software that led to different sensor outputs. The used bandwidth and output power are the same for all sensors.

Sincerely yours



Dipl.-Ing.(FH) Martin Glänzer
CTO

Difference to KY-LOC 1D.02.01:

- KY-LOC 1D.01.01
 - This will be a product with a reduced resolution. The HF signals are identical.
- KY-RAY 1D.01.02
 - This product is measuring against a passive target instead of a second device.
- KY-RAY 1D.02.02
 - This product is the same as the KY-RAY 1D.01.02 but with a higher accuracy.