

Synchroni Scientific Instruments Inc.

Date: June 14, 2024

FEDERAL COMMUNICATIONS COMMISSION
AUTHORIZATION AND EVALUATION DIVISION
7435 OAKLAND MILLS ROAD
COLUMBIA, MD 21046

We, the undersigned hereby declare, that the following radio equipment:

Applicant : Synchroni Scientific Instruments Inc.
Address : 11 Lazenby Crt, Markham ON, L6C 1S7, Canada
Product name: Wearable Bio-Sensor
Model : Synchroni Sync-Uno, Synchroni Sync-Trio, Synchroni Sync-Pento, Synchroni Sync-Octo,
Synchroni Sync-Neo
FCC ID : 2BGWZ-SYNC
FRN : 0035547132
Freq Range: 2402~2480MHz

is electrically equal to the **original equipment**:

Applicant : OYMotion Technologies Co., Ltd.
Address : Floor 6, Building 2, 222 Guangdan Road, Pudong, Shanghai PRC
Product name: Wearable Bio-Sensor
Model : Synchrony OB5000S OL, Synchrony OB5000S BL, Synchrony OB5000S SL,
Synchrony OB5000 OL, Synchrony OB5000 BL, Synchrony OB5000 SL,
Synchrony OB5000+ OL, Synchrony OB5000+ BL, Synchrony OB5000+ SL,
Synchrony OB5000Pro OL, Synchrony OB5000Pro BL, Synchrony OB5000Pro SL,
Synchrony OB5000MAX OL, Synchrony OB5000MAX BL, Synchrony OB5000MAX SL,
Synchrony OB5000S OM, Synchrony OB5000S BM, Synchrony OB5000S SM,
Synchrony OB5000 OM, Synchrony OB5000 BM, Synchrony OB5000 SM,
Synchrony OB5000+ OM, Synchrony OB5000+ BM, Synchrony OB5000+ SM,
Synchrony OB5000Pro OM, Synchrony OB5000Pro BM, Synchrony OB5000Pro SM,
Synchrony OB5000MAX OM, Synchrony OB5000MAX BM,
Synchrony OB5000MAX SM, Synchrony OB5000S OS, Synchrony OB5000S BS,
Synchrony OB5000S SS, Synchrony OB5000 OS, Synchrony OB5000 BS,
Synchrony OB5000 SS, Synchrony OB5000+ OS, Synchrony OB5000+ BS,
Synchrony OB5000+ SS, Synchrony OB5000Pro OS, Synchrony OB5000Pro BS,
Synchrony OB5000Pro SS, Synchrony OB5000MAX OS, Synchrony OB5000MAX BS,
Synchrony OB5000MAX SS

FCC ID : 2AWVV-OB5000
FRN : 0029752326
Date of Grant: 05/29/2024
Freq Range : 2402~2480MHz

Except the applicant information, manufacturer information, model and product external photos had changed, the EUT are electronically/electrically identical to the original equipment and that the original test results continue to be representative of and applicable to it.

The new product version is identical to the certified product. And the document on-file with FCC from the original certification, remain unchanged and continue to be valid.

BEST REGARDS,

Synchroni Scientific Instruments Inc.



Ziteng Ni / CEO
Name / Title