

Model difference letter

2025.03.12

Dear Reviewer,

This is a new Application from **Shenzhen brglight technology co., ltd**

Model No.:

**BRGRCR001,BRGRCR011,BRGRCR002,BRGRCR003,BRGRCR004,
BRGRCR005,BRGRCR006,BRGRCR007,BRGRCR008,BRGRCR009,
BRGRCR010,BRGRCR012,BRGRCR013,BRGRCR014,BRGRCR015,
BRGRCR016,BRGRCR017,BRGRCR018,BRGRCR019,BRGRCR020**

Use **BRGRCR001,BRGRCR011** for all tests.Their electrical circuit design, layout, components used and internal wiring are identical,Only the model name and Appearance (Number of buttons) is different. Added differential testing for radiation and conduction.

Thank you for your cooperation!

Sincerely,

Ronggui Bi

Signature

Ronggui Bi

Manager

Shenzhen brglight technology co., ltd

Room 601, Building A5,Anle Industrial Zone, No.172,Hangcheng
Avenue,Sanwei Community,Hangcheng Street, Baoan
District,Shenzhen, Guangdong