

Photographic Documentation

CUSTOMER: SIEMENS MEDICAL SOLUTIONS
EQUIPMENT: TWELVE LEAD BIOMEDICAL TELEMETRY
608-614MHZ TRANSMITTER/RECEIVER
TESTED BY: RONALD H. CROOKER
OPERATING MODE: NORMAL

DATE: 22 JULY 2002
TEST NUMBER: 1

COUPLING DEVICE: EMCO
TEST SPEC: FCC PART 15 SUBPART B AND
PART 15.242



Photograph Description: Radiated set-up

FORM CTS-PHOTO

Photographic Documentation

CUSTOMER: SIEMENS MEDICAL SOLUTIONS
EQUIPMENT: TWELVE LEAD BIOMEDICAL TELEMETRY
608-614MHZ TRANSMITTER/RECEIVER
TESTED BY: RONALD H. CROOKER
OPERATING MODE: NORMAL

DATE: 22 JULY 2002
TEST NUMBER: 1

COUPLING DEVICE: EMCO
TEST SPEC: FCC PART 15 SUBPART B AND
PART 15.242



Photograph Description: Radiated set-up

FORM CTS-PHOTO

Photographic Documentation

CUSTOMER: SIEMENS MEDICAL SOLUTIONS
EQUIPMENT: TWELVE LEAD BIOMEDICAL TELEMETRY
608-614MHZ TRANSMITTER/RECEIVER
TESTED BY: RONALD H. CROOKER
OPERATING MODE: NORMAL

DATE: 22 JULY 2002
TEST NUMBER: 2
COUPLING DEVICE: LISN
TEST SPEC: FCC PART 15 SUBPART B



Photograph Description: Conducted set-up

FORM CTS-PHOTO

Siemens Medical Solutions
Twelve Lead Biomedical Telemetry 608-614MHz Transmitter/Receiver