

1. THIS SPECIFICATION IS CONTROLLED BY CRHF PACKAGE ENGINEERING. THIS SPECIFICATION IS INTENDED TO BE THE FINAL LABEL DESIGN OUTPUT GENERATED FROM A COMBINATION OF LABEL PARTS WITHIN THE ASSOCIATED PDM BOM.
2. THESE LABEL IMAGES ARE THE COMBINATION OF THE TEMPLATE FILE, MODEL DATA RECORD, SAMPLE TRIGGER DATA AND REPRESENTATIVE VIEWS OF THE LABEL DIE LINES.
3. PRINTING REQUIREMENTS - THE LABELS ARE TO BE PRINTED USING A 300 DPI ZEBRA THERMAL TRANSFER PRINTER.
4. THE LOT OR SERIAL NUMBER, USE BY DATE, DATE OF MANUFACTURE, PIN NUMBER, ADDRESS CODE, GTIN, PRODUCT NUMBER (REF) AND BAR CODES ARE SHOWN FOR EXAMPLE ONLY. THESE FIELDS ARE VARIABLE AND ARE TO BE SUPPLIED BY MANUFACTURING AT THE TIME OF PRINTING.
5. SPACING AND TEXT WRAPPING SHOWN IN THE IMAGES MAY DIFFER SLIGHTLY FROM ACTUAL PRINTED LABEL, AS LONG AS CONTENT OF PRINTED VERSION IS ACCURATE. THE IMAGES INCLUDE REPRESENTATION OF THE PRE-PRINTED LABEL STOCK. BRANDING, PART NUMBER BAR CODES (2D OR OTHER) AND OTHER PRE-PRINTED AREAS OR ADDITIONAL STICKERS MAY CHANGE AND ARE ACCEPTABLE.
6. THE LABEL PRINTED IN MANUFACTURING MUST BE AT AN EQUAL OR LATER MODEL DATA AND TEMPLATE REVISION AS SHOWN ON THE IMAGE.
7. (FOR LABEL SET ONLY) THE TRIGGER DATA SHOWN ON THE LABEL SET IMAGE MAY NOT MATCH THE PRINTED LABEL.

SPEC TAB	MFG ADDRESS
A86279-001	Tolochenaz, Switzerland
A86279-002	Juncos, Puerto Rico, USA

DeviceImageSpec_02


Medtronic MEDTRONIC, INC. MINNEAPOLIS, MINNESOTA USA MEDTRONIC CONFIDENTIAL THIS DOCUMENT IS THE PROPERTY OF MEDTRONIC, INC. AND MUST BE ACCOUNTED FOR. INFORMATION HEREON IS CONFIDENTIAL. DO NOT REPRODUCE IT, REVEAL IT TO UNAUTHORIZED PERSONS OR SEND IT OUTSIDE MEDTRONIC WITHOUT PROPER AUTHORIZATION.	TITLE LABEL IMAGE SPEC-SOLARA QUAD CRT-P MRI SURESCAN W4TR03 US				
	LABEL IMAGE SPEC TEMPLATE: A83320 REV C				
	DATE 07-Feb-2017	SCALE N/A	ENGINEERING JUDY SALZER		
	MODEL W4TR03	Size B	PART NUMBER A86279	REV B	SHEET 1 of 3



Medtronic


Bluetooth™

Solara™ Quad CRT-P MRI SureScan™





CRT-P



SureScan™



MR Conditional



Burst+/- RA



Ramp: RA



Package contents:

1x



Implantable device



Torque wrench



Product documentation

DDD



↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate



0.3 mV

Sensitivity: RA



0.9 mV

Sensitivity: RV



→ 130 ms
→ 100 ms

A-V interval (paced and sensed)



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RA



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RV



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: LV



LV1 - LV2

Pacing polarity: LV



LV->RV
0 ms

V pacing and V-V pace delay



AUTO

Pacing polarity: RA, RV



AUTO

Sensing polarity: RA, RV



OFF

Therapies: AT, AF



AUTO

PVARP



Adaptive

Rx Only

FCFID: LFSBLEIMPLANT2

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Manufactured in: Tolochenaz, Switzerland
www.medtronic.com

STERILE EO



-18 °C
0 °F



55 °C
131 °F





M968371A001 B - M970304A001 A




Medtronic


Bluetooth™

Solara™ Quad CRT-P MRI SureScan™



CRT-P





MR Conditional

Mode **W4TR03**

SN WWW000000W

2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

REF W4TR03

Package contents:

1x

Implantable device

Torque wrench

Burst + RA

Ramp: RA

Package contents:

1x

Implantable device

Torque wrench

DDD



↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate

0.3 mV

Sensitivity: RA

0.9 mV

Sensitivity: RV

→ 130 ms
→ 100 ms

A-V interval (paced and sensed)

↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RA

↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RV

LV1 - LV2

Pacing polarity: LV

LV->RV
0 ms

V pacing and V-V pace delay

AUTO

Pacing polarity: RA, RV

Adaptive

AUTO

Sensing polarity: RA, RV

OFF

Therapies: AT, AF

AUTO

PVARP

Rx Only

FCC ID: LFSBLEIMPLANT2

STERILE

EO




Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA

Manufactured in: Tolochenaz, Switzerland

www.medtronic.com

-18 °C
0 °F



55 °C
131 °F





M968371A001 B - M970304A001 A

M968371A001 B - M970304A001 A	
	 2000-00-00 YYYY-MM-DD
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	  Solara™ Quad CRT-P MRI SureScan™  WWW000000W

Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation date:

2000-00-00

YYYY-MM-DD

REF

W4TR03

SN

WWW000000W

Medtronic

SureScan™

MR

Bluetooth®

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation date:

2000-00-00

YYYY-MM-DD

SN

WWW000000W

Medtronic

SureScan™

MR

Bluetooth®

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation date:

2000-00-00

YYYY-MM-DD

SN

WWW000000W



Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

OAE-DDDR

2x IS-1 BI **1x IS4-LLLL**



(01)00000000000000(17)000000(21)WWW000000W

REF

W4TR03

SN

 WWW000000W



2000-00-00
 YYYY-MM-DD

M968371A001 B - M970304A001 A



End Label 1



OAE-DDDR

2x

IS-1 BI

1x

IS4-LLLL

Solara™ Quad CRT-P MRI SureScan™ W4TR03

- VectorExpress™ LV Automated Test
- PhysioCurve™ Design
- OptiVol™ 2.0 Fluid Status Monitoring



REF W4TR03

SN WWW000000W



2000-00-00
YYYY-MM-DD

M968371A001 B - M970304A001 A

MR Conditional

PIN: W0W0W00000W

End Label 2

SPEC TAB
A86279-001

A86279-001


NOTE: IMAGES NOT TO SCALE
 LABEL IMAGE SPEC TEMPLATE: A83320 REV C


PART NUMBER	REV	SHEET
A86279	B	2 of 3


Medtronic



Bluetooth

Solara™ Quad CRT-P MRI SureScan™


CRT-P


SureScan


MR Conditional



RA, IS-1
RV, IS-1
LV, IS4 - LLLL

OAE-DDDR

ACCELEROMETER

Burst+: RA

Ramp: RA

Package contents:

1x

Implantable device

Torque wrench

Product documentation

Model **W4TR03**

SN **WWW000000W**

2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

REF **W4TR03**

DDD

↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate

0.3 mV

Sensitivity: RA

0.9 mV

Sensitivity: RV

→ 130 ms
→ 100 ms

A-V interval (paced and sensed)

↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RA

↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RV

↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: LV

LV1 - LV2

Pacing polarity: LV

LV->RV
0 ms

V pacing and V-V pace delay

AUTO

Pacing polarity: RA, RV

AUTO

Sensing polarity: RA, RV

OFF

Therapies: AT, AF

AUTO

PVARP

Adaptive

Rx Only

FCID: LFSBLEIMPLANT2

STERILE EO

18 °C
0 °F

55 °C
131 °F

Info

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Manufactured in: Juncos, Puerto Rico, USA
Product of USA

www.medtronic.com

M968371A001 B - M970340A001 A

Medtronic

Bluetooth™

Solara™ Quad CRT-P MRI SureScan™

CRT-P

Mode W4TR03

SN WWW000000W

2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

REF W4TR03

DDD

↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate

0.3 mV

Sensitivity: RA

0.9 mV

Sensitivity: RV

→ 130 ms
→ 100 ms

A-V interval (paced and sensed)

↑ 3.5 V
→ 0.4 ms

Amplitude and pulse width: RA

↑ 3.5 V
→ 0.4 ms

Amplitude and pulse width: RV

↑ 3.5 V
→ 0.4 ms

Amplitude and pulse width: LV

LV1 - LV2

Pacing polarity: LV

LV->RV
0 ms

V pacing and V-V pace delay

AUTO

Pacing polarity: RA, RV

Adaptive

AUTO

Sensing polarity: RA, RV

OFF

Therapies: AT, AF

AUTO

PVARP

Rx Only

FCC ID: LFSB1EIMPLANT2

STERILE

EO

55 °C
131 °F

-18 °C
0 °F

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA

Manufactured in: Juncos, Puerto Rico, USA

www.medtronic.com

M968371A001 B - M970304A001 A

M968371A001 B - M970304A001 A		
		2000-00-00 YYYY-MM-DD
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W
Medtronic W4TR03	 SureScan	 MR
		Solara™ Quad CRT-P MRI SureScan™  WWW000000W

Medtronic

Solaris™ Quad CRT-P MRI SureScan™ W4TR03

Implantation
date: 2000-00-00
YYYY-MM-DD

SureScan™

MR Conditional

REF W4TR03

SN WWW000000W

Medtronic    **Bluetooth®**

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation
date:  **2000-00-00**
SN **WWW000000W** **YYYY-MM-DD**




Medtronic    **Bluetooth®**

Solar™ Quad CRT-P MRI SureScan™ W4TR03

Implantation
date:  **2000-00-00**
YYYY-MM-DD

SN **WWW000000W**





Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

OAE-DDDR

2x IS-1 BI 1x IS4-LLLL



(01)00000000000000(17)000000(21)WWW000000W

REF W4TR03

SN WWW000000W

2000-00-00
YYYY-MM-DD

M968371A001 B - M970304A001 A

End Label 1


Medtronic

OAE-DDDR

2x **IS-1 BI** **1x** **IS4-LLLL**

Solara™ Quad CRT-P MRI SureScan™ W4TR03

- VectorExpress™ LV Automated Test
- PhysioCurve™ Design
- OptiVol™ 2.0 Fluid Status Monitoring



SureScan

MR

MR Conditional

REF W4TR03

SN WWW000000W

2000-00-00
YYYY-MM-DD

PIN: W0W0W00000W

M968371A001 B - M970304A001 A

End Label 2

SPEC TAB
A86279-002

DeviceImageSpec_02

NOTE: IMAGES NOT TO SCALE
 LABEL IMAGE SPEC TEMPLATE: A83320 REV C

REV	SHEET
B	3 of 3