

Resting ECG Analysis System
FCC ID: OU5-MAC7
IC: 4048B-MAC7

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

REVISIONS				
REV	MADE BY	DATE	ECR/ECO NO.	DESCRIPTION
2	Zhang Long	SEE MWS	SEE MWS	FCC and IC Number update.

LBL P/N: 2111601-001

Notes:

1. Use the following color callouts listed below:

- White
- Black

2. Blank label P/N: 422161-001.

UNLESS OTHERWISE SPECIFIED:
1) REMOVE ALL BURRS AND SHARP EDGES.
2) DIMENSIONS ARE IN MILLIMETERS.
X ±0.1
X.X ±0.1
X.XX ±0.01
X° ±0.5°
3) M PROFILE METRIC THREAD
TOLERANCE CLASS 6H/6g PER
ISO 9651/1.

ESSENTIAL
CHARACTERISTIC: NO

DO NOT SCALE DRAWING
CONFIDENTIAL: THIS PRINT AND THE DESIGN THEREIN REPRESENTED ARE THE PROPERTY OF GE HEALTHCARE AND MAY NOT BE USED IN ANY MANNER DETRIMENTAL TO IT'S INTEREST. THIS PRINT HAS BEEN ISSUED WITH THE UNDERSTANDING THAT IT WILL NOT BE USED FOR MANUFACTURE FOR OTHERS, AND WILL BE RETURNED UPON REQUEST.

DRAWN BY
ZHANG LONG | DATE
10-Jan-2019

ENG. APPROVED
ELECTRONICALLY APPROVED

MODEL NAME
--

 GE Healthcare

TITLE
MAC 7 WLAN LABEL


DRAWING NUMBER
2111601-001 | REV
2

SCALE 1:1 | SHEET 1 OF 1

CTQ TABLE		
QTY	FLAG	DESCRIPTION
0		CTQ FLAG: MEFI40 QUALITY PLAN REQUIRED.



INSPECTION KEY:
ITEMS INDICATED MAY REQUIRE
A SUPPLIER INSPECTION PLAN.

BALLOON KEY:  XX = FIND NUMBER, YY = QUANTITY