

Date: January 9, 2025

To whom it may concern

Cover Letter – Infinity® M540 patient monitor

IC: 5895B-MS40251

FCC ID: URA-MS40251

The host with HMN: “Infinity® M540 patient monitor” is a patient monitor with local alarm functionality. It supports wired and wireless data transmission depending on the usage.

There are the following application possibilities

- The host is locked into a holder on the patient’s bed. As soon as the monitor is locked WLAN is deactivated and data is only transmitted via LAN.
- The host can be removed from the holder by medical staff (occupational use) and can be then used as hand-held portable device. In this case WLAN is activated.

There is no possibility of wearing the device on the body.

The distance from antenna to hand is about 11.19 mm as described in the document “M540 Antenna Antenna Location.pdf”

Sincerely,



Georgios Kokovidis

Wireless Manager

Dräger Medical Systems Inc.