

CAIDENT

# iTero™ Intra-Oral Scanner

Optical Impression Device

Model: EID-2

P/N: TA15100-A

~120V/60Hz/5A

This product complies with:  
21 CFR 1040.10 & IEC 60825-1



Class 1 Laser Product

This device complies with Part 15 of FCC Rules.  
Operation is subject to the following two conditions:  
(1) this device may not cause harmful interference, and  
(2) this device must accept any interference received,  
including interference that may cause undesired operation.

S/N: 93.5-0150

Data Code: 0307

FCC ID: UZ2MA14730

Made in Israel

